



Instructions for Use

Name: Disposable Dental Box

The product Dental Box varies as follows,

- 1) Retainer Box (Large) – RB 001
- 2) Retainer Box (Small) – RB 002

Intend Use: Disposable Dental Boxes are in the form of a container to hold the denture or retainer. They are intended to be used to keep a denture or retainer clean.



Cautions:

1. Please check if the packaging is completed before use.
2. The device should be used before expiry date.

Instruction for use:

1. Put denture or retainer inside it.
2. Dispose after use to comply with regulation.

Storage: The product should be stored in dry place and keep away from acid, heat and direct sunlight during storage and transportation.

Shelf life: 5 years

Labels, Packing Logo Design:

Symbol	Introductions	Symbol	Introductions
	Lot Number / Batch Code		Do not reuse" are "single use, "Use only once
	Warnings and Precautions		non-sterile
	medical device		Manufacture Date
	Manufacturer Name & Address		Expiry Date





Sif Medical Manufacturing

	Do not use if package is damaged		Consult Instructions for use
	Mo direct sunlight		Keep dry

Manufacturer Information



Saudi Industrialists Factory for Plastic Products
Factory No.105, Zone A,
Sudair City for Industry & Business, Riyadh – Saudi Arabia
Tel: +966 112774111 / +966 553469922
Website: www.sif-manufacturing.com E-mail: info@sif-manufacturing.com

Note : Occurrence of any serious incidents should report to the manufacturer and SFDA (19999) immediately



 00966553469922

 info@sif-manufacturing.com

 BOX 88992-RIYADH 11692

 sif-manufacturing.com